

Clinical and experimental study of a new ceramic-enforced PEEK-titanium hybrid abutment placed in healed bone with an immediate loading ceramic restoration

Application of polymer-based abutments for final restorations

JOSÉ EDUARDO MATÉ-SÁNCHEZ DE VAL AND JOSÉ LUIS CALVO-GUIRADO, MURCIA



Literaturangabe

- Adell R, Eriksson B, Lekholm U, Brånemark PI, Jemt T. Long-term follow-up study of osseointegrated implants in the treatment of totally edentulous jaws. *Int J Oral Maxillofac Implants.* 1990 Winter;5(4):347-59. [Links]
- Jemt T, Lekholm U. Implant treatment in edentulous maxillae: a 5-year follow-up report on patients with different degrees of jaw resorption. *Int J Oral Maxillofac Implants.* 1995 May-Jun;10(3):303-11. [Links]
- Rocci A, Rocci M, Rocci C, Scoccia A, Gargari M, Martignoni M, Gottlow J, Sennerby L. Immediate loading of Brånemark system TiUnite and machined-surface implants in the posterior mandible, part II: a randomized open-ended 9-year follow-up clinical trial. *Int J Oral Maxillofac Implants.* 2013 May-Jun;28(3):891-5. [Links]
- Cibirka RM, Nelson SK, Lang BR, Rueggeberg FA. Examination of the implant-abutment interface after fatigue testing. *J Prosthet Dent.* 2001 Mar;85(3):268-75. [Links]
- Goodacre CJ, Bernal G, Rungcharassaeng K, Kan JY. Clinical complications with implants and implant prostheses. *J Prosthet Dent.* 2003 Aug;90(2):121-32. [Links]
- Ribeiro CG, Maia MLC, Scherrer SS, Cardoso AC, Wiskott HWA. Resistance of three implant-abutment interfaces to fatigue testing. *J Appl Oral Sci.* 2011 Aug;19(4):413-20. [Links]
- Strub JR, Gerds T. Fracture strength and failure mode of five different single-tooth implant-abutment combinations. *Int J Prosthodont.* 2003 Mar-Apr;16(2):167-71. [Links]
- Tripodakis AP, Strub JR, Kappert HF, Witkowski S. Strength and mode of failure of single implant all-ceramic abutment restorations under static load. *Int J Prosthodont.* 1995 May-Jun;8(3):265-274. [Links]
- An in vitro pilot study of abutment stability during loading in new and fatigue-loaded conical dental implants using synchrotron-based radiography. Rack T, Zabler S, Rack A, Riesemeier H, Nelson K. *Int J Oral Maxillofac Implants.* 2013 Jan-Feb;28(1):44-50. doi: 10.11607/jomi.2748. Erratum in: *Int J Oral Maxillofac Implants.* 2013 Jul-Aug;28(4):1002. Rack, Tanja [corrected to Rack, Tatjana].
- Assessment of human gingival fibroblast interaction with dental implant abutment materials. Rutkunas V, Bukelskiene V, Sabaliuskas V, Balciunas E, Malinauskas M, Baltriukiene D. *J Mater Sci Mater Med.* 2015 Apr;26(4):169. doi: 10.1007/s10856-015-5481-8. Epub 2015 Mar 25.
- Steinberg EL, Rath E, Schlaifer A, Chechik O, Maman E, Salai M. Carbon fiber reinforced PEEK Optima - A composite material biomechanical properties and wear/debris characteristics of CF-PEEK composites for orthopedic trauma implants [Internet]. *J Mech Behav Biomed Mater* [cited 2014 07 4]. 2013 Jun;17:221-8. Available from: <http://www.sciencedirect.com/science/article/pii/S1751616112002573>. DOI: 10.1016/j.jmbbm.2012.09.013. [Links]
- Invisio Biomaterial Solutions. PEEK-OPTIMA POLYMER. Product manual [homepage on the Internet]. 2004 [cited 2012 Feb 15]. Available from: <http://www.invisio.com>. [Links]
- Kurtz SM, Devine JN. PEEK biomaterials in trauma, orthopedic, and spinal implants. *Biomaterials* 2007;28:4845-4869.
- Toth JM, Wang M, Estes BT, Scifert JL, Seim HB 3rd, Turner AS. Polyetheretherketone as a biomaterial for spinal applications. *Bio- materials* 2006;27:324-334.
- Tetelman ED, Babbush CA. A new transitional abutment for immediate aesthetics and function. *Implant Dent* 2008;17:51-58.
- Santing HJ, Meijer HJ, Raghoobar GM, Ozcan M. Fracture strength and failure mode of maxillary implant-supported provisional single crowns: A comparison of composite resin crowns fabricated directly over PEEK abutments and solid titanium abutments. *Clin Implant Dent Relat Res* 2012;14:882-889.
- Bayer S, Komor N, Kramer A, Albrecht D, Mericske-Stern R, Enkling N. Retention force of plastic clips on implant bars: A randomized controlled trial. *Clin Oral Implants Res* 2012;23:1377-1384.
- Tannous F, Steiner M, Shahin R, Kern M. Retentive forces and fatigue resistance of thermoplastic resin clasps. *Dent Mater* 2012; 28:273-278.
- Kurtz SM, Devine JN. PEEK biomaterials in trauma, orthopedic, and spinal implants. *Biomaterials* 2007;28:4845-4869.
- Toth JM, Wang M, Estes BT, Scifert JL, Seim HB 3rd, Turner AS. Polyetheretherketone as a biomaterial for spinal applications. *Bio- materials* 2006;27:324-334.
- Tetelman ED, Babbush CA. A new transitional abutment for immediate aesthetics and function. *Implant Dent* 2008;17:51-58.
- Hallmann L, Mehl A, Sereno N, Haemmerle CH. The improvement of adhesive properties of PEEK through pre-treatments. *Appl Surf Sci* 2012;258:7213-7218.